


Ball Attachments

The titanium **Nova** Ball Attachment provides firm retention and stabilization of the overdenture when it is used with tissue-supported implant-retained overdentures.

The titanium **Nova** Ball Attachment is compatible with internal hex.

Standard

Ø 4.00 mm



Ø 4.00 mm

Cat. No.	BAT-0010	BAT-0020	BAT-0030	BAT-0040	BAT-0050	BAT-0060
Platform Height	1.00 mm	2.00 mm	3.00 mm	4.00 mm	5.00 mm	6.00 mm
Length	8 mm	9 mm	10 mm	11 mm	12 mm	13 mm
Platform	3.75 mm	3.75 mm	3.75 mm	3.75 mm	3.75 mm	3.75 mm

Covers and Caps

Cat. No.	HUS - 0000	SIL 0001	SIS 0002	SIH 0003
				
Placement	Ball Cover	Cap	Cap	Cap
Material	Metal	Silicone	Silicone	Silicone
Finish		Standard	Soft	Hard