


Standard

Ø 4.76 mm




Ø 4.76 mm

Cat. No.	ZRX-0910	ZRX-1110	ZRX-1111
Height	1 mm	2 mm	3 mm
Length	9 mm	10 mm	11 mm
Platform	3.75 mm	3.75 mm	3.75 mm

Nova manufactures top quality prosthetic abutments. Use with the SRR 9.4 mm standard screw and the standard a 1.25mm hex driver.

Angulated 15°

Ø 4.76 mm




Ø 4.76 mm

Cat. No.	ZRX-1115-10	ZRX-1215-20	ZRX-ZRX-1315-30
Angle α	15°	15°	15°
Height	1 mm	2 mm	3 mm
Length	9 mm	10 mm	11 mm
Platform	3.75 mm	3.75 mm	3.75 mm

Angulated 25°

Ø 4.76 mm



Ø 4.76 mm

Cat. No.	ZRX-1125-10	ZRX-1125-20
Angle α	25°	25°
Height	1 mm	2 mm
Length	9 mm	10 mm
Platform	3.75 mm	3.75 mm

Nova manufactures top quality prosthetic abutments. Use with the SRR 9.4 mm standard screw and the standard a 1.25mm hex driver.