



Nova Implants are tapered internal hex implants, designed to enable easy insertion while supporting excellent initial stability. The variable thread design enables self tapping, thus providing superb solutions for a variety of bone conditions.

All our implants are made of the titanium alloy Ti-6Al-4V ELI, in accordance with ASTM-F136-02. Titanium has the outstanding capability to integrate almost completely with the bone. This unique characteristic coupled with its mechanical advantages of strength & endurance makes titanium the ideal implant material. In addition it can be precisely machined to ensure a range of implants that meets the requirements of modern state of the art Prosthetics. All our implants can be ordered with the regular SLA surface treatment or the new **RBM** more advanced surface treatment.

Nova Implants multi-design features offers an excellent solution for immediate placement and immediate loading. One platform fits all implant diameters. Packaging has been designed for quick identification and easy opening.

Nova Implants adhere to the highest international standards of manufacturing and have the ISO standards Institute and CE approval.

Clinical Benefits:

- Better control during placement.
- Greater initial placement stability.
- Self drilling, super self-tapping feature.
- Self-condensing.
- Self-directing in line of insertion.

Product Identification

NAME	DIAMETER	LENGTH
PSI -	3.75 -	L8

Color Identification - implants



Ø 3.3 mm



Ø 3.3 mm

Cat. No.	PSI-3.3-L8	PSI-3.3-L10	PSI-3.3-L11.5	PSI-3.3-L13	PSI-3.3-L16
Length	8 mm	10 mm	11.5 mm	13 mm	16 mm
Surface Treatment	SLA / RBM	SLA / RBM	SLA / RBM	SLA / RBM	SLA / RBM

Ø 3.75 mm



Ø 3.75 mm

Cat. No.	PSI-3.75-L8	PSI-3.75-L10	PSI-3.75-L11.5	PSI-3.75-L13	PSI-3.75-L16
Length	8 mm	10 mm	11.5 mm	13 mm	16 mm
Surface Treatment	SLA / RBM	SLA / RBM	SLA / RBM	SLA / RBM	SLA / RBM

Ø 4.2 mm



Ø 4.2 mm

Cat. No.	PSI-4.2-L8	PSI-4.2-L10	PSI-4.2-L11.5	PSI-4.2-L13	PSI-4.2-L16
Length	8 mm	10 mm	11.5 mm	13 mm	16 mm
Surface Treatment	SLA / RBM	SLA / RBM	SLA / RBM	SLA / RBM	SLA / RBM

Ø 5.0 mm



Ø 5.0 mm

Cat. No.	PSI-5-L8	PSI-5-L10	PSI-5-L11.5	PSI-5-L13	PSI-5.0-L16
Length	8 mm	10 mm	11.5 mm	13 mm	16 mm
Surface Treatment	SLA / RBM	SLA / RBM	SLA / RBM	SLA / RBM	SLA / RBM

Ø 6.0 mm



Ø 6.0 mm

Cat. No.	PSI-6-L8	PSI-6-L10	PSI-6-L11.5	PSI-6-L13
Length	8 mm	10 mm	11.5 mm	13 mm
Surface Treatment	SLA / RBM	SLA / RBM	SLA / RBM	SLA / RBM

Color Identification - drills



Implant Procedures

